

2025-11-26
EURAMET Project 1782

**Comparison in terms of reference air kerma rate for HDR Ir-192 using
well-type chambers**

Identifier in the BIPM key comparison database (BIPM KCDB):
EURAMET.RI(I)-K8.x

Technical Protocol
(approved by EURAMET TC-IR and CCRI(I))

Pilot laboratory: STUK (reporting) and SSM (measurements)

Contact persons:

Reetta Nylund

STUK – Radiation and Nuclear Safety Authority

Jokiniemenkuja 1

01370 Vantaa, Finland

Tel. +358 40 152 0941

e-mail: reetta.nylund@stuk.fi

Linda Persson

Swedish Radiation Safety Authority (SSM)

Riksmätplatsen

Solna Strandväg 122

17154 Solna, Sweden

Tel: +46 70 917 66 71

e-mail: linda.persson@ssm.se

Contents

1	Object and participants	1
1.1	Objective of the comparison	1
1.2	Participants.....	1
2	Transfer instruments	2
3	Measurement Procedure	4
3.1	Radiation qualities and quantities	4
3.2	Reference conditions	4
3.3	¹⁹² Ir source	4
3.4	Reference value.....	5
3.5	Determination of the calibration coefficient.....	5
3.6	Uncertainty budgets.....	6
3.7	Reporting the results	6
3.8	Evaluation of the results.....	7
4	Course of comparison	8
4.1	Transport and time schedule	8
5	References	10
	Appendix I: Complete addresses of the participants used for shipments	12
	Appendix II. Pictures of transfer chambers.....	14

1 Objective and participants

The national metrology laboratories for dosimetry quantities of the Netherlands (VSL, the National Metrology Institute of the Netherlands), Sweden (SSM, Swedish Radiation Safety Authority), Finland (STUK, Radiation and Nuclear Safety Authority), Norway (DSA, Norwegian Radiation and Nuclear Safety Authority), Italy (ENEA, Agenzia Nazionale per le nuove tecnologie), Serbia (VINS, Vinča Institute of Nuclear Sciences - National Institute of the Republic of Serbia, University of Belgrade), Greece (EEAE, Greek Atomic Energy Commission), Croatia (RBI, Ruđer Bošković Institute) and the International Atomic Energy Agency (IAEA) have agreed to perform a comparison in terms of reference air kerma rate (RAKR) for HDR ^{192}Ir using well-type chambers. All participating laboratories, apart from VSL, are secondary standards dosimetry laboratories (SSDLs) for the respective quantity. IAEA is the only SSDL to maintain its own ^{192}Ir source, while the other SSDLs perform calibrations at hospitals.

In the project two well-type transfer chambers will be circulated among participants, and each laboratory will report calibration coefficients and their expanded uncertainties for these chambers in terms of reference air kerma rate. VSL as a primary standards dosimetry laboratory will provide a link to the BIPM.RI(I)-K8 comparison throughout this comparison.

This technical protocol prepared by the laboratories specifies the procedure to be followed in this dosimetry comparison. The technical protocol is prepared according to the BIPM technical protocol for the ongoing BIPM key comparison BIPM.RI(I)-K8. The purpose of a comparison is to compare the calibration capabilities of the participating laboratories, not to require each participant to adopt precisely the same conditions of measurement. The protocol, therefore, specifies the procedures necessary for the comparison, e.g. reference conditions, but not the procedures used in the calibration of the laboratories being compared.

1.1 Objective of the comparison

The objective of the comparison is to support the ionising radiation CMCs of SSM, IAEA, EEAE, RBI, VINS, ENEA, STUK, and DSA, in the dosimetry branch for the quantity reference air kerma rate (RAKR) for HDR ^{192}Ir .

1.2 Participants

Table 1 lists the participants. In Appendix I, the complete contact details for the participants are presented. Table 2 presents the traceability of participating laboratories in terms of RAKR for ^{192}Ir .

Table 1. Participants of the project in the order of measurements.

Institute	Country	Contact person	e-mail
SSM (pilot, measurements)	Sweden	Linda Persson	linda.persson@ssm.se
VSL (reference value)	Netherlands	Grazia Brazzorotto Jacco de Pooter	gbrazzorotto@vsl.nl jdpooter@vsl.nl

IAEA	International, Austria	Zakithi Msimang	z.msimang@iaea.org
EEAE	Greece	Argiro Boziari	argiro.boziari@eeae.gr
RBI	Croatia	Luka Bakrac	lbakrac@irb.hr
VINS	Serbia	Milos Zivanovic	milosz@vin.bg.ac.rs
ENEA	Italy	Alessia Ciccotelli	alessia.ciccotelli@enea.it
STUK (pilot, reporting)	Finland	Reetta Nylund	reetta.nylund@stuk.fi
DSA	Norway	Per Otto Hetland	per.otto.hetland@dsa.no

Table 2. Traceability of calibrations at the participating laboratories in terms of reference air kerma rate for HDR ^{192}Ir .

Institute	^{192}Ir Traceability	Type of standard
SSM	PTB	Secondary
VSL	VSL	Primary
IAEA	PTB	Secondary
EEAE	LNE-LNHB	Secondary
RBI	PTB via IAEA	Secondary
VINS	PTB via IAEA	Secondary
ENEA	VSL	Secondary
STUK	NPL	Secondary
DSA	PTB via IAEA	Secondary

2

Transfer instruments

Two well-type ionization chambers of the same type and manufacturer will be used as transfer instruments for this comparison. Two chambers will be used to minimize the risk of potential transfer instrument breakage during the comparison. If a chamber suffers a failure during the comparison, the comparison will be continued with the other chamber, listed in this protocol. If there are results for a particular chamber from at least a half of the participating laboratories, results will be reported in the final publication. The chambers are owned by different laboratories participating in the comparison and these laboratories have stability data for these chambers available, which can be used if needed to evaluate the chamber stability during the comparison. No electrometer is circulated, and laboratories shall use their own electrometers and cables for performing the measurements during the comparison and ensure traceability for their current measurement and high-voltage setting applied to the ion chamber. The details of the measurement equipment (e.g. afterloader used) shall be specified when reporting results.

The technical details of the chambers are listed in table 3.

Table 3. Technical data of the transfer chambers (chamber characteristics according to IAEA TRS492 and manufacturer's information).

Chamber type (owner of chamber)	Standard Imaging HDR 1000 Plus (SSM)	Standard Imaging HDR 1000 Plus (IAEA)
Serial number	A062722	#A953382
Geometry	well-type	well-type
Active volume	245 cm ³	245 cm ³
Height	15.6 cm	15.6 cm
Diameter	10.2 cm	10.2 cm
Insert height	12.1 cm	12.1 cm
Insert diameter	3.5 cm	3.5 cm
Accompanying insert	Standard Imaging HDR Ir source holder 70010, labelled Br2 RMP06	Standard Imaging HDR Ir source holder 70010
Polarising voltage of a chamber	+300 V on collector (central) electrode, 0 V on chamber wall (collecting negative charge) (if +300 V on collector is not available: -300 V on chamber wall, 0 V on collector electrode)	+300 V on collector (central) electrode, 0 V on chamber wall (collecting negative charge) (if +300 V on collector is not available: -300 V on chamber wall, 0 V on collector electrode)
Connector type	BNC (BNT 2-lug)	BNC (BNT 2-lug)
Nominal source holder reference points	Black dot on the insert and dot on the top of the chamber housing	Black dot on the insert and dot on the top of the chamber housing
Sweet spot	Shall be determined separately by each participant (nominal maximum response at about 50 mm from the bottom of the chamber insert)	Shall be determined separately by each participant (nominal maximum response at about 50 mm from the bottom of the chamber insert)
Other remarks	Figure A1	Figure A2

3 Measurement procedure

The dosimetry laboratories are expected to ensure that their reference standard is in perfect working order prior to the comparison. When the participant receives the transfer chambers, they shall perform a visual check for any damage and verify correct functioning prior to any additional measurements. If it seems that the chambers are broken, the participant should contact the piloting laboratories STUK and SSM to discuss further actions.

Each participant will proceed following their own calibration procedure(s) according to their quality management system to determine the calibration coefficients of the transfer chambers in terms of reference air kerma rate. Whenever possible, all calibration results shall be given corrected to Elekta Flexisource HDR Ir-192 and plastic catheter, which is the reference for this comparison. This can be performed by applying the source model correction factor e.g. according to the method in IAEA TRS-492 and the data in TRS-492 Table 8, and the participant needs to report how the conversion is performed. An uncorrected value can be given if the conversion to Flexisource is not possible. Then the participant is responsible for considering the effect of the estimated correction in the uncertainty budget. Furthermore, each laboratory may add any other correction factors needed for calibrations and these correction factors shall be reported when reporting results.

3.1 Radiation qualities and quantities

The radiation quality used in the comparison is ^{192}Ir . The quantity used for the comparison is reference air kerma rate defined according to ICRU85a. Primary standards used in the comparison and as source of traceability (Table 2) should have the ICRU90 recommendations regarding the reference data for primary standards implemented.

3.2 Reference conditions

The calibration coefficients for the transfer chambers should be given in terms of reference air kerma rate in units of Gy/C, determined at the point of maximum chamber response (sweet spot) and corrected to standard conditions of air temperature and pressure; $T_0 = 293.15 \text{ K}$ and $P_0 = 101.325 \text{ kPa}$. Corrections for relative humidity (RH = 50%) will not be applied for the measurements performed between 20% and 80% RH. Each laboratory will use their own equipment to measure environmental conditions and ensure traceability for those measurements.

3.3 ^{192}Ir source

Each participant may choose their usual source construction/type and activity. As the well chamber responses are dependent on the ^{192}Ir source design, the source reference code, manufacturer and activity should be reported. Furthermore, the results should be corrected into same geometry as the reference value for this comparison using the source model correction factor.

3.4 Reference value

VSL as linking laboratory to BIPM.RI(I)-K8 will provide the link to the BIPM.RI(I).K8 KCRV (Key Comparison Reference Value) for this comparison. VSL will proceed according to the same procedures as for the latest comparisons.

VSL will perform reference measurements using an Elekta Flexitron afterloader and plastic catheter for Standard Imaging chambers.

The reference values will be calculated separately for each chamber and the values reported by SSDLs will be compared to the reference value by calculating a ratio between calibration coefficients of the SSDL and linking laboratory (VSL), e.g.

$$K_{RAKR,SSDL}/K_{RAKR,LINK}.$$

3.5 Determination of the calibration coefficient

Each laboratory details their own procedure or refers to international practices/guidance followed when performing the calibration. Several of SSDLs will perform measurements at hospitals and will follow their own procedures for this.

Typically, the laboratory measures the position of maximum chamber response (sweet spot) by stepping a brachytherapy source through the well chamber. This position shall be reported in the comparison report, for instance as the distance from the tip of the (plastic) catheter (which is pushed to the bottom of the borehole in the chamber holder) and the centre of the Ir-192 source at the sweet spot, including the length of the catheter.

Typically, the SSDL establishes a reference air kerma rate of an ^{192}Ir HDR source in accordance with their own procedure following an equation such as:

$$\dot{K}_{RAKR} = N_{\dot{K}_{RAKR,SSDL}} I_{SSDL} \quad (1)$$

where $N_{\dot{K}_{RAKR,SSDL}}$ is the calibration coefficient for reference air kerma rate determined by a primary standards dosimetry laboratory (PSDL) for a specific HDR ^{192}Ir source model, and where I_{SSDL} is the ionisation current measured by the SSDL with an electrometer system with traceability to electrical standards. In accordance with TRS-492, I_{SSDL} is corrected to reference conditions of air temperature and pressure ($T=20\text{ °C}$ and $P=101.325\text{ kPa}$) (k_{Tp}), and if needed for relative humidity, chosen for the comparison. Furthermore, the relevant source model correction factor, k_{sm,sm_0} , should be applied to $N_{\dot{K}_{RAKR,PSDL}}$ if the type of source used at the SSDL is different to the source used for the calibration of the well-type chamber at the PSDL. For the corrections to I_{SSDL} , a laboratory shall proceed according to their own procedure and may include e.g. the electrometer correction factor, correction for leakage, correction for distance, correction for volume, correction for ion recombination etc. All corrections used shall be reported in addition to the final results and they all contribute to the uncertainty budget.

Each SSDL calibrates the transfer chamber for reference conditions using their usual calibration procedure to derive the calibration coefficient for the transfer chamber, $N_{RAKR,tr}$, which is defined as:

$$N_{RAKR,tr} = \frac{\dot{K}_{RAKR}}{I_{tr}} \quad (2)$$

where $\dot{K}_{RAKR}$ is the reference air kerma rate from equation (1), and where I_{tr} is the signal from the transfer chamber measured by the SSDL at the specific reference polarity stated in table 3 using their own electrometer systems with traceability to electrical standards. Like I_{SSDL} , I_{tr} is corrected to standard conditions of air temperature and pressure, and if needed for relative humidity, chosen for the comparison. Similar corrections as for I_{SSDL} should be applied to I_{tr} .

VSL will determine calibration coefficients in the same way as the SSDLs, except that realization of the quantity RAKR is performed at VSL using their primary method to realize reference air kerma, as described in the respective comparison reports (Kessler *C. et al.*, 2024 and van Dijk *et al.*, 2004).

3.6 Uncertainty budgets

In addition to calibration coefficients each participant shall provide a detailed measurement uncertainty budget considering all influence parameter and correction factors for each calibration coefficient. Each participant shall describe the main components of the uncertainty including correction factors in the budget in the level of one standard uncertainty and provide the final expanded combined uncertainty, $k=2$. The detailed uncertainty budget shall be provided in accordance with the Guide to the Expression of Uncertainties in measurements (JCGM, 2008) with corresponding confidence level. Components of the uncertainty budget shall be provided as relative values [%]. The Excel report sheet includes an example form for the uncertainty budget, into which each laboratory is recommended to add components according to their procedures.

3.7 Reporting the results

The pilot laboratory (STUK) will send their results to the CCRI Executive Secretary Vincent Gressier (vincent.gressier@bipm.org) within 6 weeks of completing their measurements. Other participants will send their results (calibration coefficients and uncertainty budgets) with information of the measuring set-up to STUK (reetta.nylund@stuk.fi) after STUK has sent their measurements to the CCRI Executive Secretary. As the STUK will be one of the last participants performing the measurements, all other participants will keep their results confidential until then. SSM will deliver data for chamber stability within 2 weeks of completing their stability measurements, after which STUK will begin data analysis. If a participant has not sent their results (calibration coefficient and uncertainty budgets) by the due date, the laboratory will be excluded from the comparison.

A common Excel template for reporting the results will be provided to each participant in addition to this technical protocol.

Before the draft A is delivered to participants, the pilot laboratory (STUK) will confirm with all participants that they will participate using the given results. If there is not enough information available, e.g., uncertainty budgets don't include all needed components to estimate/calculate degrees of equivalence for the comparison, the pilot

laboratory (STUK) reserves the right to contact the participant to obtain the details. In this case a participant is expected to answer within two weeks to the pilot laboratory to keep the comparison on track.

3.8 Evaluation of the results

After the reporting pilot laboratory (STUK) has received all results, the results of the SSDLs will be evaluated in comparison with VSL's results, as VSL provides a link to the KCRV of this comparison. The results will be analysed for each transfer chamber. Degrees of equivalence will be calculated in terms of comparison results (calibration coefficients and uncertainties) according to CCRI(I)/17-09 instructions. In general, the degree of equivalence of each SSDL, with respect to the KCRV, is evaluated as follows separately for each transfer chamber:

$$R_{SSDL} = \frac{N_{RAKR,SSDL}}{N_{RAKR,VSL}} R_{VSL,BIPM} \quad (3)$$

in which $R_{VSL,BIPM}$ represents the results of VSL in the BIPM.RI(I)-K8 key comparison (Kessler C et al, 2024).

The variance of R_{SSDL} is:

$$u_{R,SSDL}^2 = (u_{SSDL}^2 + u_{BIPM}^2 - \sum_j f_j^2 (u_{SSDL,j}^2 + u_{BIPM,j}^2)) + u_{tr}^2 + u_{LINK}^2 \quad (4)$$

in which f_j are weighting factors related to correlating components, which will be evaluated during the analysis.

In equation (4) u_{tr} combines the stability of the transfer chambers over the period of the comparison and the variation in the ratios for a specific chamber. u_{tr} will be calculated based on standard uncertainties of R_{SSDL} (from equation (5)) for each laboratory according to CCRI(I)/17-09. Additionally, the pilot laboratory (SSM) will perform several measurements for each transfer chamber and u_{tr} may also be calculated and adjusted based on these (e.g. if SSM's repeatability data is significantly higher than u_{tr} , i.e. more than two standard deviations). If needed, stability data from chamber owners may be used to evaluate chamber stability. If u_{tr} is significantly higher (more than two standard deviations) for a specific chamber, that chamber will be excluded from the comparison.

In equation (4) u_{VSL} represents the uncertainty of non-statistical components, which are not cancelling out via the linking mechanism. To ease estimation of these, VSL's measurement conditions in this comparison are as close as possible to those used in the VSL-BIPM comparison. However, a possible uncertainty due to intra-type variation of sources might be added into the linking mechanism.

For the case of the 2 transfer chambers which will be circulated in this comparison, equation (3) becomes:

$$R_{SSDL} = \left(\frac{N_{RAKR,VSL}}{N_{RAKR,BIPM}} \right) \frac{1}{2} \sum_{j=1}^2 \frac{N_{RAKR,SSDL,j}}{N_{RAKR,VSL,j}} \quad (5)$$

Equation (5) will eventually be used for the analysis, if both transfer chambers are used for the analysis of results.

The degree of equivalence for the SSDL is:

$$D_{SSDL} = R_{SSDL} - 1 \quad (6)$$

and its expanded uncertainty is $U_{SSDL}=2u_{R,SSDL}$.

Further details for data analysis may be discussed among the participants on the basis of the Draft A report.

In the reporting of the results, document “CIPM MRA-G-11: Measurement comparisons in the CIPM MRA, Guidelines for organizing, participating and reporting” will be followed.

4 Course of comparison

4.1 Transport and time schedule

The laboratory should make all the arrangements for safe transport of the transfer standards once measurements have been completed. Each participating institute is responsible for its own costs regarding the measurements and transportation of the transfer instruments to the next institute. The standards won't be insured by the pilot, and each participant is responsible for the good care of the chambers in their facilities and good packing of the chambers for the subsequent shipment. Each participant is responsible for the chambers within their country and onward shipment until the receiving laboratory has received the equipment. Shipment shall be made using a courier. It is recommended to take photographs of the chambers before the shipment.

The laboratories outside of the EU (Serbia, Norway) will perform their measurements in the neighbouring EU-country as onsite calibration (as they would perform calibrations at a hospital in their home countries) and thus, the transfer chambers won't be shipped outside of the EU. Laboratories outside of the EU are responsible for their own travel to neighbouring countries as well as transport of their equipment.

IAEA will ship to SSM at IAEA cost at the beginning of the comparison for stability measurements and SSM will provide shipping of the chambers back to IAEA at the end of the comparison.

The transfer standards are packed in a Pelican protection box together with a complete description of the devices (i.e. technical protocol) including information about the manufacturer, type, serial number, size, weight and technical data needed for their operation. The package is 55 cm x 26 cm x 48 cm and it weighs 15.7 kg.

Table 4 summarises the proposed schedule of the comparison measurements and table 5 summarises the course of the comparison events. Comparison measurements will begin before the final approval of the technical protocol due to the availability of ¹⁹²Ir.

Table 4. Prosed schedule for the comparison measurements. The measurements have been started before the full acceptance of the technical protocol.

Institute	Measuring period	Date of chamber leaving to next participant
SSM, stability check	October 2025	October 31st, 2025
VSL	November 2025	Nov 30 th , 2025
IAEA	December 2025	Dec 31st, 2025
EEAE	January 2026	Jan 31st, 2026
RBI	February 2026	Feb 28th, 2026
VINS	February 2026	
ENEA	March 2026	Mar 30 th , 2026
STUK	April 2026	Apr 30 th , 2026
DSA	May 2026	
SSM, measurements and stability check	May-June 2026	June 15 th , 2026

Table 5. Proposed comparison events.

Event	Due date
IAEA ships a transfer standard to SSM	October 3 rd , 2025
SSM ships both transfer standards to VSL	October 31 st , 2025
Comparison measurements (see table 4)	November 2025
STUK sends results to BIPM	May 31 st , 2026
Participants shall not send their results before to STUK!	
Final due date for results to STUK	June 30 th , 2026
Chamber stability measurement reported (from SSM to STUK)	June 30 th , 2026
Draft A delivered to participants	September 15 th , 2026
Comments by the participants on draft A	October 10 th , 2026

Draft B available	December 1 st , 2026
Final report available	Depending upon comments by evaluators

5

References

BIPM, BIPM.RI(I)-K8 Technical Protocol, version 6.0, July 2024

[https://www.bipm.org/kcdb/comparison/doc/download/285/bipm.ri\(i\)-k8_technical_protocol%20.pdf](https://www.bipm.org/kcdb/comparison/doc/download/285/bipm.ri(i)-k8_technical_protocol%20.pdf)

CCRI(I)/17-09, Burns DT, Butler D, Updated report on the evaluation of degrees of equivalence in regional dosimetry comparisons, May 2017

CIPM MRA-G-11, Measurement comparisons in the CIPM MRA

Guidelines for organizing, participating and reporting, version 1.1, 18/01/2021

<https://www.bipm.org/documents/20126/43742162/CIPM-MRA-G-11.pdf/9fe6fb9a-500c-9995-2911-342f8126226c>

IAEA, International Atomic Energy Agency. Dosimetry in Brachytherapy – An International Code of Practice for Secondary Standards Dosimetry Laboratories and Hospitals. Technical Reports Series No. 492, June 2023.

https://www-pub.iaea.org/MTCD/Publications/PDF/DOC-010-492_web.pdf

ICRU, ICRU85a (Revised) Fundamental Quantities and Units for Ionizing Radiation, <https://www.icru.org/report/fundamental-quantities-and-units-for-ionizing-radiation-icru-report-85a-revised/>

ICRU, ICRU Report 90, Key Data For Ionizing-Radiation Dosimetry: Measurement Standards And Applications, <https://www.icru.org/report/icru-report-90-key-data-for-ionizing-radiation-dosimetry-measurement-standards-and-applications/>

JCGM (Joint Committee for Guides in Metrology). Evaluation of measurement data – Guide to the expression of uncertainty in measurement. JCGM 100:2008, GUM 1995 with minor corrections. First edition, September 2008.

Kessler C, Jansen B J, de Pooter JA, Aviles Lucas P, Key comparison BIPM.RI(I)-K8 of high dose-rate ¹⁹²Ir brachytherapy standards for reference air kerma rate of the VSL and the BIPM, *Metrologia* **61** 06010, 2024

DOI 10.1088/0026-1394/61/1A/06010

<https://iopscience.iop.org/article/10.1088/0026-1394/61/1A/06010>

van Dijk E, Kolkman-Deurloo IK, Damen, PMG, Determination of the reference air kerma rate for ¹⁹²Ir brachytherapy sources and the related uncertainty, *Med. Phys.*, **31**, 2826-2833, 2004

Appendix I: Complete addresses of the participants used for shipments

SSM / Sweden

Swedish Radiation Safety Authority (SSM)
 Riksmätplatsen
 Solna Strandväg 122
 171 54 SOLNA
 Sweden
 Contact person: Linda Persson
 Tel: +46 70 917 66 71
 e-mail: Linda.Persson@ssm.se

VSL/the Netherlands

VSL, National Metrology Institute
 Thijsseweg 11
 2629 JA DELFT
 The Netherlands
 Contact person: Grazia Brazzorotto/ Jacco de Pooter
 Tel: +31615356332 / +31631119896
 e-mail: gbrazzorotto@vsl.nl / jdpooter@vsl.nl

IAEA/Austria

IAEA Laboratories Seibersdorf
 Dosimetry Laboratory
 Friedensstrasse 1
 A-2444 Seibersdorf
 Austria
 Contact persons: Zakithi Msimang and Joao Cardoso
 Tel: +43-1-2600-21660 / 28328
 e-mail: z.msimang@iaea.org / j.cardoso@iaea.org

EEAE/Greece

Greek Atomic Energy Commission, Ionizing Radiation Calibration Laboratory
 Patriarxou Grigoriou and Neapoleos
 15341 Agia Paraskevi
 Attika
 Greece
 Contact person: Argiro Boziari
 Tel: +302106506765
 e-mail: argiro.boziari@eeae.gr

RBI / Croatia

Ruder Bošković Institute
Secondary standard dosimetry laboratory
Bijenička cesta 54
10 000 Zagreb
Croatia
Contact person: Luka Bakrač
Tel: +385 99 6989 169
e-mail: lbakrac@irb.hr

ENEA / Italy

ENEA-INMRI Istituto Nazionale di Metrologia delle Radiazioni Ionizzanti
C.R. Casaccia
Via Anguillarese, 301,
00123 Roma RM
Italy
Contact person: Alessia Ciccotelli
Tel: +39 0630483357
e-mail: alessia.ciccotelli@enea.it

STUK / Finland

Radiation and Nuclear Safety Authority (STUK)
Radiation Metrology Laboratory
Jokiniemenkuja 1
01370 Vantaa
Finland
Contact person: Reetta Nylund
Tel: +358 401520941
e-mail: reetta.nylund@stuk.fi

Appendix II. Pictures of transfer chambers.

Standard Imaging HDR 1000 Plus A062722




Figure A1. Standard Imaging HDR 1000 Plus, s/n A062722. The picture on the right shows the correct alignment of the Standard Imaging HDR Iridium Source Holder, model 70010, labelled Br2 RMP06.

Standard Imaging HDR 1000 Plus A953382



Figure A2. Standard Imaging HDR 1000 Plus, s/n A953382. The other transfer chamber is exactly similar as another chamber (the source holder in figure is not correct holder).