

13.10.2016



High-pressure qNMR spectroscopy in condensed- and gas-phase towards determination of impurities and compositions of gas mixtures

Klas Meyer, Heinrich Kipphardt, Michael Maiwald

Online NMR Spectroscopy as Reference

Hyphenation to reactors and devices

500 MHz
Online NMR



reactors

Complex reacting multicomponent mixtures

- quantitative information
- qualitative („chemical“) information
- highest possible dispersion

Direct coupling to reactors

- fast and direct sample transfer
- automation
- extended pressure range




Process NIR



Process UV/VIS




Process Raman



Process NMR


How does NMR work ?

NMR stands for **N**uclear **M**agnetic **R**esonance



Radio Frequency pulse

For $I=1/2$




- e.g. ^1H , ^{13}C
- High dispersion of signals
- Absolute comparable method
- Quantitative information by integration

Quantitative NMR spectroscopy

Purity assessment of raw materials

- *iso*-Pentane for production of Primary Reference Gas Mixtures




- Orientational results from ^{13}C -NMR
- Comparison with estimated GC-MS

Component	$x_{\text{NMR}} / \%$	$x_{\text{GC-MS}} / \%$
<i>iso</i> -Pentan	99.794	99.771
<i>n</i> -Pentane	0.170	0.187
<i>neo</i> -Pentane	0.036	0.038

Quantitative NMR spectroscopy

Purity assessment of raw materials

- *neo*-Pentane: problem of available purity limited to 2.0




- Results from 7 spectra
- Uncertainty ~ 120 ppm (n/n)

Component	$x_{\text{NMR}} / \%$	$u(x_{\text{NMR}}, k=2) / \%$
<i>neo</i> -Pentane	98.945	0.012
Cyclobutane	1.055	0.012

High-Pressure NMR spectroscopy

Application in gas and condensed phase

Gas




- Higher density of nuclei
- Direct connection to gas cylinder
- Pressure ≤ 20 MPa (up to 100 MPa possible)

→ Increasing sensitivity,
 → ^{13}C -NMR in the gas phase



Liquid




- Sampling from piston cylinders
 - Circulation for homogenization
- Direct quantitation in liquids
 → No sample preparation



Gas-phase NMR spectroscopy


Primary Reference Gas Mixtures

- Methane-based mixture with components up to C₅ (BAM-G412, H1-8K)
- Available filling pressure: ~3 MPa



¹H-NMR

32 Scans, time ~ 10 min



¹³C-NMR

128 Scans, time ~ 68 min


512 Scans, time ~ 273 min

→ Several hours of measurement time is needed, especially for ¹³C-NMR


Gas-phase NMR spectroscopy

Long-term measurements

- No deuterated substances for field-frequency-lock
- Signal drifts during experiment
- Standard procedure of accumulation not suitable
- Alignment required before accumulation of spectra



without alignment




after alignment

Gas-phase NMR spectroscopy

Accumulation procedure

- Accumulation of single spectra in frequency domain
- Preprocessing required → Phasing, Baseline correction
- Accumulated spectrum can be treated like standard NMR data




Gas-phase NMR spectroscopy

Primary Reference Mixture BAM-G412

Degree of Equivalence

$$\frac{x_{NMR}}{x_{cert}} - 1$$


- Linefit with gaussian/lorentzian functions for signal area
- Relative quantification with 100%-approach
- Expanded uncertainty ($k = 2$):
 - Contribution of repeated measurement
 - Contribution of repeated data analysis
- Comparison with certified gravimetric value



	CH ₄	C ₂ H ₆	C ₃ H ₈	<i>n</i> -C ₄ H ₁₀	<i>iso</i> -C ₄ H ₁₀	CO ₂
$x(\text{certified})$	97.283	1.028	0.253	0.208	0.204	0.965
$u_x(\text{CRM cat.}, k = 2)$	±0.0482	±0.0040	±0.0020	±0.0016	±0.0016	±0.0045
$t = 47.2 \text{ h}$	97.354	1.051	0.198	0.200	0.203	0.994
$t = 133.1 \text{ h}$	97.401	0.978	0.222	0.217	0.210	0.972
$t = 161.6 \text{ h}$	97.241	1.028	0.210	0.208	0.216	1.097
$x(\text{NMR})$	97.332	1.019	0.210	0.208	0.210	1.021
$u_x(\text{combined}, k = 2)$	±0.0951	±0.0432	±0.0141	±0.0098	±0.0071	±0.0772
Degree of Equivalence	-8.877E-05	-0.0096	-0,1702	0.0031	0.0265	0.0576

Gas-phase NMR spectroscopy

Degree of Equivalence BAM-G412




- Certified value mostly within uncertainty range of NMR measurements
- Unresolved deviation in case of propane

NMR spectroscopy in condensed phase

Idea and Motivation



- High-purity hydrocarbon mixtures commonly provided in piston cylinders
- In use for key comparisons, e.g., EURAMET.QM-S6/1195 (NPL)



- Avoiding of evaporation for homogeneous withdrawal of content
- Enrichment of light volatile components in NMR cell expected

NMR spectroscopy in condensed phase

Experimental setup for circulation





- Brazed concentric lines 1/32" in 1/8"
- Circulation of low viscosity fluids like liquefied gases
- PFA tubing in active region of spectrometer
- Currently tested up to 6 MPa (higher expected)

NMR spectroscopy in condensed phase

Results on liquefied hydrocarbons

- Commercial mixture, linear alkanes C₃ to C₇ (~20 mol-% each)




- Enrichment and Homogenization observed
- Unexpected high scatter in repeatability

NMR spectroscopy in condensed phase

Results on liquefied hydrocarbons




- Comparison of spectra pretreatment and evaluation methods
- Uncertainties ($k=2$) ~600–2000 ppm estimated based on repeatability



NMR spectroscopy in condensed phase

Results on liquefied hydrocarbons

- Two sample cylinders from CCQM-K119 key comparison




- All signals can be used for integration
- Better repeatability in static mixture

NMR spectroscopy in condensed phase

Results on liquefied hydrocarbons



- Four samplings from cylinder with 7-10 spectra on each sample
- Uncertainties between 900 and 2800 ppm (n/n) → problem 100%-method



NMR spectroscopy in condensed phase

Construction of a piston cell

- Displacement element for increasing pressure of gas samples
- Measurements in expanded liquid phase and supercritical fluids
- Sapphire window for visual observation and optical methods




- Pressure up to 60 MPa, temperature planned up to 130 °C
- Prototype available

NMR spectroscopy in condensed phase

Construction of a piston cell

- Mounting in a reinforced frame with guide rails
- Spindle drive unit for piston movement, $x = 0.25 \text{ mm / turn}$



- Long-term pressure tests at 60 MPa (water)
- Tests with Helium showed good pressure stability

Summary

Purity assessment with qNMR

- Determination of impurities in raw materials for PRG production
- Identification of compounds, no prior calibration necessary

NMR in the gas phase

- High-pressure setup for gas sampling up to 20 MPa
- Long-term ^{13}C -NMR measurements on primary reference gas mixtures
- Development of accumulation procedure with drift correction

NMR in condensed phase

- Setup for studies of liquefied hydrocarbons from piston cylinders
- Investigations on test mixture and samples of CCQM-K119 comparison

Acknowledgement




Paul Brewer
Lucy Culleton

Thank you for your attention !


Gas-phase NMR spectroscopy

Advantages of high pressure

- mixture of Propene 3.5 and Methane 5.5
- $T = 27\text{ }^\circ\text{C}$, $p = 10.6\text{ MPa}$



$^1\text{H-NMR}$
1 Scan \cong 5 s



$^{13}\text{C-NMR}$
8 Scans \cong 750 s

→ Significant savings in measurement time with increase of pressure