

Joint Research Centre

the European Commission's in-house science service

*Serving society
Stimulating innovation
Supporting legislation*


JCTLM Member Activities: IRMM

H. Schimmel, G. Auclair, I. Zegers
E. Monogioudi, J. Kuhlmann

SID Unit

Sèvres, November 30th 2015

The JRC inside the European Commission



The Joint Research Centre (JRC)



European Commission's in-house science service

Established in 1957 (Euratom Treaty)


Supporting EU policies with independent, evidence-based scientific and technical support

7 Institutes (Directorates) in 6 locations

~ 3.000 staff (2015)

Budget €393 million annually, plus €73 earned income

EU Legislative procedure



Revision and implementation of IVD Directive (98/79/EC) / Regulation

- *Proposal published 2012, IRMM consulted.*
- *Since then under scrutiny and amended by the Council of the European Union, European Parliament and European Commission (information see http://ec.europa.eu/growth/sectors/medical-devices/regulatory-framework/revision/index_en.htm).*
- *In 2015 several 'trilogue' negotiations between Council of the European Union, European Parliament and European Commission to sort out remaining issues (progress being made).*
- *IRMM represented in IVD technical group, consisting of competent authorities representatives addressing technical issues concerning implementation of the IVD Directive \ Regulation*


Overview activities

- *CRM development concentrating on:*
 - *CRMs for enzymatic catalytic activity (renewals)*
 - *HbA1c calibration solutions (changed format)*
 - *Serum proteins (B2M)*
 - *CRMs for Alzheimer's disease diagnostics (Abeta42)*
 - *CRMs for autoimmune disease diagnostics*
- *Provide advice to IFCC and its committees and working groups*
- *Provide input to ISO TC212 (revision ISO 17511)*

CRMs for HbA1c

- *At the final stage*
- *Gravimetric mixtures of pure HbA1c and HbA0*
- *6 levels at approximately 0, 30, 60, 90, 120 and 150 mmol/mol*
- *Can be directly used for the calibration of the IFCC reference measurement procedure for HbA1c*

CRM for B2M



ERM-DA470k/IFCC

Released in 2008 and certified for 12 proteins, freeze-dried


β -2-microglobulin

marker for multiple myeloma and lymphoma


2.17 ± 0.07 mg/L (k=2)

CRM for B2M

ERM-DA470k/IFCC using kit calibrators



ERM-DA470k/IFCC using calibrant provided by IRMM





CRMs for autoimmune diagnostics (IFCC WG-HAT)


- *Two layer standardisation problem:*
 - *heterogeneous response between individuals to variable antigens (at least antigen to be defined)*
 - *standardisation of the quantitative measurement of the IgG fraction binding to a defined antigen*

CRMs for autoimmune diagnostics – correlation studies


Manufacturers	Platforms & Kits
AESKU	AESKULISA® MPO
Bio-Rad	Anti-MPO EIA
	Bioplex 2200
Eurodiagnostica	MPO ANCA DIASTAT®
	Capture MPO-ANCA Wieslab®
	MPO-ANCA Wieslab®
Euroimmun	Anti-MPO ELISA
	Myeloperoxidase (MPO) (pANCA) (IgG)
IMMCO	Anti-MPO ELISA
Inova	QUANTA Lite MPO SC ELISA
	MPO BIOFLASH
	QUANTA Lite® MPO IgG
Menarini	Menarini
Orgentec	ORG 519 Anti-MPO (pANCA
Phadia	Varelisa™ MPO ANCA
	EliA MPO ^S
Orgentec	Anti-MPO ORG 519



CRMs for autoimmune diagnostics - commutability studies



CRMs for autoimmune diagnostics – value assignment




CRMs for autoimmune diagnostics – value assignment to calibrant





ERM-DA470k/IFCC

**Nephelometry
Turbidimetry**




**Purified
anti-MPO IgG**

CRMs for autoimmune diagnostics – value assignment to calibrant



$$[\text{anti-MPO IgG}] = \frac{\text{slope}_{\text{anti-MPO IgG}}}{\text{slope}_{\text{ERM-DA470k/IFCC}}} \times [\text{ERM-DA470k/IFCC}]$$

CRMs for autoimmune diagnostics – value assignment to matrix CRM



**Purified
anti-MPO IgG**

**anti-MPO IgG
methods**




ERM-DA476/IFCC

CRMs for autoimmune diagnostics – value assignment to matrix CRM

Linear function (IRMM- Wieslab CAP MPO IU)


Average Values



■ anti-MPO IgG ERM-DA476/IFCC

Logarithmic function IMMCO

Average values



■ anti-MPO IgG ERM-DA476/IFCC

$$[\text{ERM-DA476/IFCC}] = \frac{\text{slope}_{\text{ERM-DA476/IFCC}}}{\text{slope}_{\text{anti-MPO IgG}}} \times [\text{anti-MPO IgG}]$$

CRMs for autoimmune diagnostics

Method	Method values		Common calibrant	
	Value	Unit	Value	Unit
Wieslab® Capture MPO-ANCA	40	IU/mL	92	mg/L
ImmuLisa™	356	IU/mL	68	mg/L
Wieslab® Capture MPO-ANCA	45	IU/mL	88	mg/L
Wieslab® MPO-ANCA	76	IU/mL	73	mg/L
EliA MPO ^S	77	IU/mL	93	mg/L
Autoimmune EIA Anti-Myeloperoxidase	78	U/mL	99	mg/L
AESKULISA MPO	205	U/mL	91	mg/L
Anti-Myeloperoxidase ELISA (IgG)	>200	RU/mL	75	mg/L
QUANTA Flash MPO	98	CU	78	mg/L

CRMs for autoimmune diagnostics



Certified value [mg/L]	U_{CRM} ($k = 2$) [mg/L]
84	9

- **Will the reference materials solve all problems?** **NO**
 - outliers, different selectivity, calibration strategies
- **What can the materials contribute to?** **An anchor point**
 - reduction of batch to batch variation
 - long-term stability
 - better equivalence of measurement results


CRMs for autoimmune diagnostics – IRMM and IFCC WG-HAT

	anti-MPO	anti-PR3	anti-B2GP
Raw materials	✓	✓	✓
Correlation	✓	✓	✓
Commutability	✓	✓	✓
IgG purified	✓	✓	✓
Processing	✓	✓	✓
Homogeneity	✓		✓
Stability	✓		✓
Value assignment	✓		On-going

✓ -done at IRMM


✓ -done by other partners of the working group (Silvia Pierangeli, Rohan Willis)

CRMs for Alzheimer's diagnostics – correlation and commutability




Bjerke et al, CCLM, in print


CRMs for Alzheimer's diagnostics – material processing



Drift during filling (aCSF)



problem at beginning and end of filling (aCSF)




Results from last trial (CSF)

$S_{\text{between vial}} = 1.9 \%$
(between vial variability)

by optimizing processing conditions

CRMs for Alzheimer's diagnostics – traceability chain



Acknowledgements:

- *Many colleagues at IRMM who contributed (project responsables, laboratory and processing staff)*
- *Members of the IFCC Working Groups and Committees IRMM is collaborating with*
- *Collaboration partners (researchers and IVD manufacturers)*